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Key Benefits

- Provides 3D imaging in the interventional suite to enhance decision making and guidance
- Supports accurate assessment of vascular pathologies by providing high-resolution 3D reconstructions of small vessels and lesions
- Enhances understanding of vascular anatomy for interventional treatment planning and procedural outcome verification.

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Enhancing 3D functionality

Visualizing the complex spatial relationship between critical and branching vessels often involves several sequential 2D (DSA) acquisitions and radiation dose for the patient. SmartCT Angio offers a 3D-RA (3D Rotational Angiography) acquisition protocol that provides extensive 3D visualization of anatomy and vessels based on a single contrast-enhanced rotational angiogram. Its high-resolution 3D reconstructions provide critical information about depth and the relationship of one vessel to another to support the accurate assessment of anatomy and vasculature.

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With SmartCT Angio, complex anatomy such as aneurysms, complex anatomy, or tortuous vessel structures can be assessed in three dimensions. This enhances the chances of delineating the neck of aneurysms, for example, and its shape and relationship to adjacent arteries. It also enhances the assessment of complex congenital heart disease anatomy and its relationship to adjacent structures.

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Combined with the unique whole body coverage of the X-ray system, specifically designed for 3D imaging, SmartCT Angio can cover cerebral, abdominal, cardiac, and peripheral vasculature as well as other anatomy.

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4 step Guidance.

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1. Room setup

2. Proper 3D protocol with corresponding suggestion of injection protocol (when applicable)

3. Collision free Zero dose table iso-centering

4. When to press and release the acquisition button

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Image Acquisition

Image acquisition is performed with the Rotational Angiography feature of the X-ray system with the flexibility to position the C-arm in either head or side (not F12) position.

C-arm in head position: scan range of 240 degrees with a rotation speed up to 55 degrees/sec.

C-arm in side position: scan range of 180 degrees with a rotation speed up to 30 degrees/sec.

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3D Vessel Reconstruction

The rotational run is automatically transferred and displayed as a 3D vessel model: with the Real-Time digital link (option) 125 images are reconstructed into a 3-dimensional model within seconds. Additional reconstructions, using the Reconstructive Zooming Technique, can be performed as well.

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Workflow

Step by step acquisition guidance

Automated 3D-RA process from 3D acquisition to 3D Viewing,
3D at touch screen module,

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Calibration

3D-RA calibrations are performed by Philips Customer Support.

3D-RA calibration data are stable over at least 6 months' time.

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Viewing

Real Time user interface.

Philips' CRM (Contrast Resolution Management) Technology.

Image rendering:

Volume/Surface Rendering,

MIP,

Average

Gradient rendering,

MPR (Multi-Planar Reformatting),

unlimited distance measurements calculated in the same volume, including "Quick measurement".

Volume calculation

Lesion segmentation,

Annotation,

Reconstructive Zooming Technique,

Subtraction of reconstructed volumes,

Set grey values WW/WL,

Store/Recall of user defined projections.

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Transfer to:

Optional Hard Copy unit (DICOM Print),

DICOM compatible device, supported are DICOM XA, DICOM SC, DICOM CT and DICOM 3D,

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SmartCT Vessel Analysis allows easy inspection of the vessel and device positioning with straightened, curved and cross-section reformats to support treatment planning. The curved MPR view allows you to see the whole vessel segment on one plane. The straightened reformat view of the vessel segment, where the curvature is extracted from the vessel, while preserving the longitudinal and angular position, contains a graph showing the vessel diameter along the segment. The straightened cross-section view displays an indication of the minimum and maximum diameters at the pointer location as you move it over the curved, reformat or straightened reformat view. You can choose your preferred rendering to enhance visibility of guidewires and the stretched vessel view allows you to measure the diameter of the vessel/lumen and the length of the segment/stenosis at three locations. Ring landmarks can be used to mark feeder vessels to aid navigation.

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Revealing hidden structures

The complexity of interventional procedures lies in the fact that every person's pathology is unique. Visualization in three dimensions is therefore vital to aid decision making by the clinician. Rotational angiography provides real-time 3D impressions of complex vasculature and the coronary artery tree. Rotational Angio can be used to quickly determine the projection angle for treatment.

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Specifications

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Rotational Angio acquires multiple projections with just one contrast injection via a fast rotational scan of the region of interest. A rotational scan is possible both with the X-ray systems in the side position (ceiling mounted systems) and in the head position, providing the flexibility to perform procedures virtually from head to toe.

C-arm in side position:

Max. rotation Speed: 30 degrees/s

Max. rotation Angle: 180 degrees

C-arm in head position:

Max. rotation Speed: 55 degrees/s

Max. rotation Angle: 240 degrees

Max. Frame speeds are given by the frame speed specifications of the system configuration.

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The very high movement speed allows using less contrast, whereas the very wide rotation range provides a complete evaluation of the anatomy.

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The stand is designed for a very high mechanical stability. It offers precise positioning and high reproducibility, assuring you of high quality images and excellent subtraction studies. Rotational Angio results are available on the X-ray system.

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Operation of Rotational Angiography is straight forward: the procedure is selected, set up and executed virtually in a matter of seconds, supporting high patient throughput.

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A set of dedicated acquisition programs is available on the touch screen module and can be selected at the touch of a button. The Rotational Angio is controlled from the exposure hand- or footswitch

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